


BOWMAN LACHRYMAL PROBE 2/3

Réf. IS04023



BOWMAN LACHRYMAL PROBE 2/3



Réf. IS04023

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030142
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,2	Width packaging (cm):	
Length (cm) :	12,0	Length packaging (cm) :	
Height (cm) :	0,2	Height packaging (cm) :	
Weight (kg) :	0,004	Weight packaging (kg) :	0,004
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:05 - This article is available on our website : www.gsh-med.fr with the following reference IS04023

All characteristics are given for information only and may be changed