

WIRE CUTTER FORCEPS DIA.1,5MM IN 16CM T.C.

Réf. IP80423



WIRE CUTTER FORCEPS DIA.1,5MM IN 16CM T.C.



Réf. IP80423

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040044132
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:29 - This article is available on our website : www.gsh-med.fr with the following reference IP80423

All characteristics are given for information only and may be changed