

MULTI-ACTION PLIER/CUTTER 47 CM DIAM MAX 6 MM

Réf. IP80847



MULTI-ACTION PLIER/CUTTER 47 CM DIAM MAX 6 MM



Réf. IP80847

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040127989
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,0	Width packaging (cm):	12,0
Length (cm) :	47,2	Length packaging (cm) :	47,2
Height (cm) :	2,0	Height packaging (cm) :	2,0
Weight (kg) :	0,426	Weight packaging (kg) :	0,426
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:29 - This article is available on our website : www.gsh-med.fr with the following reference IP80847

All characteristics are given for information only and may be changed