


## AMALGAM PLUGGERS SMOOTH END, FOC NO.0

Réf. IDFAFOC0



# AMALGAM PLUGGERS SMOOTH END, FOC NO.0



Réf. IDFAFOC0

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe Ir               |               |
| Marquage                   | CE        | Barcode                 | 3661040015576 |
| Guarantee (month):         | 24        | Customs Code            | 9018.9090     |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 0,5       | Width packaging (cm):   | 0,5           |
| Length (cm) :              | 17,6      | Length packaging (cm) : | 17,6          |
| Height (cm) :              | 0,5       | Height packaging (cm) : | 0,5           |
| Weight (kg) :              | 0,022     | Weight packaging (kg) : | 0,022         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



30/03/2022 15:34 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IDFAFOC0