

AMALGAM PLUGGERS SMOOTH END, FOC NO.1

Réf. IDFAFOC1



AMALGAM PLUGGERS SMOOTH END, FOC NO.1



Réf. IDFAFOC1

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040015583
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,7	Width packaging (cm):	0,7
Length (cm) :	17,8	Length packaging (cm) :	17,8
Height (cm) :	0,5	Height packaging (cm) :	0,5
Weight (kg) :	0,024	Weight packaging (kg) :	0,024
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:34 - This article is available on our website : www.gsh-med.fr with the following reference IDFAFOC1

All characteristics are given for information only and may be changed