

## AMALGAM PLUGGERS SMOOTH END, FOC NO.2

Réf. IDFAFOC2



# AMALGAM PLUGGERS SMOOTH END, FOC NO.2



Réf. IDFAFOC2

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040015590
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,0	Width packaging (cm):	1,0
Length (cm) :	17,7	Length packaging (cm) :	17,7
Height (cm) :	0,7	Height packaging (cm) :	0,7
Weight (kg) :	0,024	Weight packaging (kg) :	0,024
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:34 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IDFAFOC2

All characteristics are given for information only and may be changed