

SINGLE-ENDED AMALGAM CARRIER DIA. 2.5 MM

Réf. IDPASI25



SINGLE-ENDED AMALGAM CARRIER DIA. 2.5 MM



Réf. IDPASI25

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040015811
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,0	Width packaging (cm):	2,0
Length (cm) :	15,5	Length packaging (cm) :	15,5
Height (cm) :	1,0	Height packaging (cm) :	1,0
Weight (kg) :	0,051	Weight packaging (kg) :	0,051
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:36 - This article is available on our website : www.gsh-med.fr with the following reference IDPASI25