

STAGNARA GOUGE (NON SLIPPING HANDLE) 30 CM X 10 MM

Réf. IG09410



STAGNARA GOUGE (NON SLIPPING HANDLE) 30 CM X 10 MM



Réf. IG09410

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040018997
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,5	Width packaging (cm):	8,5
Length (cm) :	35,7	Length packaging (cm) :	35,7
Height (cm) :	1,5	Height packaging (cm) :	1,5
Weight (kg) :	0,113	Weight packaging (kg) :	0,113
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference IG09410

All characteristics are given for information only and may be changed