

ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 1 LFLOVAC DISPOSABLE - CA-MI

Réf. BCFLOV304

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE **II**



ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 1 LFLOVAC DISPOSABLE - CA-MI



Réf. BCFLOV304

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040197234
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	19,0	Height packaging (cm) :	
Weight (kg) :	2,700	Weight packaging (kg) :	2,700
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



02/06/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference BCFLOV304

ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 1 LFLOVAC DISPOSABLE - CA-MI



Réf. BCFLOV304

COMPOSITION

Code	Descriptive	Quantity
BCASPI300	Askir 118 Electrical Suction Unit, 1L Jar in Makrolon CA-MI	1
BCASPI102	Container Flovac 1L (reusable) for Aspiret, Askir 20 and 30	1
BCASPI99	DISPOSABLE FLOVAC® LINER, 1000 ML for Aspiret, Askir 20, 30	1
BCASPI135	Set of tubes Ø8x14mm+connectors Flovac ASKIR 36BR/C30/118/NEWHOSPIVAC	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI suction units	1