

## ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 2L FLOVAC DISPOSABLE - CA-MI

Réf. BCFLOV305

CA-MI

GARANTIE  
2  
ANS  
YEARS  
GUARANTEE

CE  
0123

DISPOSITIF MEDICAL  
MEDICAL DEVICE  
CLASSE II



# ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 2L FLOVAC DISPOSABLE - CA-MI

Réf. **BCFLOV305**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040197241
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	19,0	Height packaging (cm) :	
Weight (kg) :	2,700	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



31/03/2022 00:19 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCFLOV305

# ASKIR 118 ELECTRICAL SUCTION UNIT, JAR 2L FLOVAC DISPOSABLE - CA-MI



Réf. BCFLOV305

## COMPOSITION

Code	Descriptive	Quantity
BCASPI300	Askir 118 Electrical Suction Unit, 1L Jar in Makrolon CA-MI	1
BCASPI103	Container Flovac 2L Aspiret, Askir 20, 30, C30BR 31841	1
BCASPI100	Disposable Flovac liner 2L for Aspiret, Askir 20, 30, C30BR	1
BCASPI135	Set of tubes Ø8x14mm+connectors Flovac ASKIR 36BR/C30/118/NEWHOSPIVAC	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI suction units	1