

FLOVAC SET 1L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET

Réf. BCASPI500

CA-MI



GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE **II**



FLOVAC SET 1L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET



Réf. **BCASPI500**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040191676
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



31/03/2022 00:15 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI500

FLOVAC SET 1L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET



Réf. BCASPI500

COMPOSITION

Code	Descriptive	Quantity
BCASPI102	Container Flovac 1L (reusable) for Aspiret, Askir 20 and 30	1
BCASPI99	DISPOSABLE FLOVAC® LINER, 1000 ML for Aspiret, Askir 20, 30	1
BCASPI111	Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI suction units	1