

FLOVAC SET 2L FOR ASKIR 118 ET 36BR



Réf. BCASPI503



FLOVAC SET 2L FOR ASKIR 118 ET 36BR

CA-MI

Réf. BCASPI503

FEATURES			
Medical device	No		
Marquage		Barcode	3661040191706
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,844
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85

31/03/2022 00:15 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI503

All characteristics are given for information only and may be changed

2/3

FLOVAC SET 2L FOR ASKIR 118 ET 36BR**Réf. BCASPI503****COMPOSITION**

Code	Descriptive	Quantity
BCASPI103	Container Flovac 2L Aspiret, Askir 20, 30, C30BR 31841	1
BCASPI100	Disposable Flovac liner 2L for Aspiret, Askir 20, 30, C30BR	1
BCASPI135	Set of tubes Ø8x14mm+connectors Flovac ASKIR 36BR/C30/118/NEWHOSPIVAC	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI suction units	1