

NEW ASPIRET ELECTRICAL SUCTION UNIT, 2L AUTOCLAVABLE POLICARBONATE JA

Réf. BCPOLY02



CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II
B

NEW ASPIRET ELECTRICAL SUCTION UNIT, 2L AUTOCLAVABLE POLICARBONATE JA

CA-MI

Réf. BCPOLY02

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe IIa | |
| Marquage | CE0123 | Barcode | 3661040194745 |
| Guarantee (month): | 24 | Customs Code | 90189084 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 23,0 | Width packaging (cm): | |
| Length (cm) : | 19,0 | Length packaging (cm) : | |
| Height (cm) : | 16,0 | Height packaging (cm) : | |
| Weight (kg) : | 2,200 | Weight packaging (kg) : | |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -20 | Temperature max | 45 |
| Humidity min | 0 | Humidity max | 85 |



31/03/2022 00:17 - This article is available on our website : www.gsh-med.fr with the following reference BCPOLY02

All characteristics are given for information only and may be changed

2/3

NEW ASPIRET ELECTRICAL SUCTION UNIT, 2L AUTOCLAVABLE POLICARBONATE JA



Réf. BCPOLY02

COMPOSITION

| Code | Descriptive | Quantity |
|----------|---|----------|
| BCASPI01 | New Aspiret electrical suction unit 1L autoclavable polycarbonate jar | 1 |
| BCASPI81 | Jars Poli 2L for Aspiret, Askir 20,30 and C30BR | 1 |