

## MUCOSITY ASPIRATOR CA-MINEW ASKIR 20 BOTTLE / BAG UU FLOVAC 1L

Réf. BCFLOV10

CA-MI

GARANTIE  
2  
ANS  
YEARS  
GUARANTEE

CE  
0123

DISPOSITIF MEDICAL  
MEDICAL DEVICE  
CLASSE  
II



# MUCOSITY ASPIRATOR CA-MINEW ASKIR 20 BOTTLE / BAG UU FLOVAC 1L



Réf. **BCFLOV10**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040200019
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	18,0	Height packaging (cm) :	
Weight (kg) :	2,500	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	1	Humidity max	85



31/03/2022 00:17 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCFLOV10

All characteristics are given for information only and may be changed

# MUCOSITY ASPIRATOR CA-MINEW ASKIR 20 BOTTLE / BAG UU FLOVAC 1L



Réf. BCFLOV10

## COMPOSITION

Code	Descriptive	Quantity
BCASPI06	New Askir 20 electrical suction unit, 1L autoclavable polycarbonate j	1
BCASPI102	Container Flovac 1L (reusable) for Aspiret, Askir 20 and 30	1
BCASPI99	DISPOSABLE FLOVAC® LINER, 1000 ML for Aspiret, Askir 20, 30	1
BCASPI111	Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI succion units	1