

MUCOSITY ASPIRATOR CA-MI NEW ASKIR 30 BOTTLE / BAG UU FLOVAC 2L

Réf. BCFLOV17

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEED

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE II



MUCOSITY ASPIRATOR CA-MI NEW ASKIR 30 BOTTLE / BAG UU FLOVAC 2L

Réf. BCFLOV17

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040200026
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	18,0	Height packaging (cm) :	
Weight (kg) :	3,600	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



31/03/2022 00:18 - This article is available on our website : www.gsh-med.fr with the following reference BCFLOV17

MUCOSITY ASPIRATOR CA-MI NEW ASKIR 30 BOTTLE / BAG UU FLOVAC 2L



Réf. BCFLOV17

COMPOSITION

Code	Descriptive	Quantity
BCASPI12	CA-MI New Askir 30 Ekectric suction unit polycarbonate jars 1L	1
BCASPI103	Container Flovac 2L Aspiret, Askir 20, 30, C30BR 31841	1
BCASPI100	Disposable Flovac liner 2L for Aspiret, Askir 20, 30, C30BR	1
BCASPI111	Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI succion units	1