

MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 2L

Réf. BCFLOV30

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II



MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 2L

Réf. **BCFLOV30**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040197227
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	23,0	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	16,0	Height packaging (cm) :	
Weight (kg) :	2,200	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



31/03/2022 00:28 - This article is available on our website : www.gsh-med.fr with the following reference BCFLOV30

MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 2L



Réf. BCFLOV30

COMPOSITION

Code	Descriptive	Quantity
BCASPI01	New Aspiret electrical suction unit 1L autoclavable polycarbonate jar	1
BCASPI103	Container Flovac 2L Aspiret, Askir 20, 30, C30BR 31841	1
BCASPI100	Disposable Flovac liner 2L for Aspiret, Askir 20, 30, C30BR	1
BCASPI111	Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30	1
BCASPI105	Round spacer for Askir 20 and 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI504	IFU for CA-MI succion units	1