

## PULSE OXIMETER NONIN PALMSAT 2500 + CHILD SENSOR WITHOUT ALARM CLOCK

Réf. DNOXYM17

 **NONIN**



**CE**  
0123



# PULSE OXIMETER NONIN PALMSAT 2500 + CHILD SENSOR WITHOUT ALARM CLOCK



## Réf. DNOXYM17

PalmSAT oximeter without alarm. Supplied with a pediatric articulated sensor.

Palmsat pulse oximeters provide accurate monitoring of oxygen saturation and heart rate. Compact, precise, easy to use, they will meet all your monitoring needs: for emergency care, operating room, intensive care and home care.

Nonin Medical's PureSAT pulse oximetry technology uses an intelligent filter that provides accurate results even under movement, low vascularization or difficult saturation conditions.

Digital red LED display with 3 digits and 7 segments.

72-hour data storage memory.

SpO2 measurement scale: 0 to 100%.

Pulse measurement scale: 18 to 300 BPM. accuracy  $\pm 3\%$ .

Temperature Operates from  $-20^{\circ}$  to  $+50^{\circ}\text{C}$ .

Humidity Operates from 10 to 90%.

Power supply: 4 x 1.5 V AA alkaline batteries (Palsat autonomy 2500 80 hours, 2500A 60 hours) or rechargeable NiMH battery charger base (approximately 40 hours of operation) sold separately ref. 39524.

Dimension: W 7 x H 13.8 x D 3.2 cm

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040208824
Guarantee (month):	36	Customs Code	9018909010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:22 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DNOXYM17

All characteristics are given for information only and may be changed