

VACUUM CLEANER CA-MI NEW HOSPIVAC 350, AP FULL 5, 2X5 LITERS AUTOCLAV

Réf. BCASPI507

CA-MI

GARANTIE
3
ANS
YEARS
GARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE **II**



VACUUM CLEANER CA-MI NEW HOSPIVAC 350, AP FULL 5, 2X5 LITERS AUTOCLAV



Réf. BCASPI507

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661041009949
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:05 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI507