

POWER ADAPTER CABLE FOR ASKIR 36BR

CA-MI

Réf. BCASPI173



POWER ADAPTER CABLE FOR ASKIR 36BR



Réf. BCASPI173

FEATURES			
Medical device	No		
Marquage		Barcode	3661040197210
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,010
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



15/11/2022 12:45 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI173