


PAIR OF TRAINING ELECTRODS FOR DEFIB. FRED PA-1 SCHILLER


Réf. SSACCE008

SCHILLER

The Art of Diagnostics



PAIR OF TRAINING ELECTRODS FOR DEFIB. FRED PA-1 SCHILLER



Réf. **SSACCE008**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661041009208
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	19,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :	0,020	Weight packaging (kg) :	0,020
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:15 - This article is available on our website : www.gsh-med.fr with the following reference SSACCE008