

INSTALLATION FOR AIVIA S/100/METAL CUPBOARD

Réf. SSACCE005


SCHILLER

The Art of Diagnostics



INSTALLATION FOR AIVIA S/100/METAL CUPBOARD

Réf. **SSACCE005**



FEATURES			
Medical device	No		
Marquage		Barcode	3661041009178
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:15 - This article is available on our website : www.gsh-med.fr with the following reference SSACCE005