


AC/DC ADAPTOR FOR AIVIA


Réf. SSACCE013

SCHILLER

The Art of Diagnostics



AC/DC ADAPTOR FOR AIVIA



Réf. **SSACCE013**

FEATURES			
Medical device	No		
Marquage		Barcode	3661041009253
Guarantee (month):	3	Customs Code	901889010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	7,5
Length (cm) :		Length packaging (cm) :	10,5
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,265
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:16 - This article is available on our website : www.gsh-med.fr with the following reference SSACCE013