

TRANSPORT BAG FOR FRED PA-1 SCHILLER-STANDARD QUALITY

Réf. SSACCE019


SCHILLER

The Art of Diagnostics



TRANSPORT BAG FOR FRED PA-1 SCHILLER-STANDARD QUALITY

Réf. **SSACCE019**



FEATURES			
Medical device	No		
Marquage	Barcode	3661041009284	
Guarantee (month):	3	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	27,0
Length (cm) :		Length packaging (cm) :	27,0
Height (cm) :		Height packaging (cm) :	11,0
Weight (kg) :	0,500	Weight packaging (kg) :	0,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:16 - This article is available on our website : www.gsh-med.fr with the following reference SSACCE019

All characteristics are given for information only and may be changed