

## PROBE COVER PC7200 RADIANT, BOX OF 200 PCS

Réf. CRAPC7200



CE

DISPOSITIF MEDICAL  
MEDICAL DEVICE  
CLASSE I

# PROBE COVER PC7200 RADIANT, BOX OF 200 PCS

Réf. **CRAPC7200**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	4710810451546
Guarantee (month):	0	Customs Code	9025900090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,0
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	3,5
Weight (kg) :		Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	50
Humidity min		Humidity max	85



31/03/2022 03:19 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference CRAPC7200