

EASY 3 BLACK ADULT SPHYGMOMANOMETER WITH 3 VELCRO CUFFS (AD-CHILD-XL)

Réf. DS99325_N



EASY 3 BLACK ADULT SPHYGMOMANOMETER WITH 3 VELCRO CUFFS (AD-CHILD-XL)



Réf. DS99325_N

FEATURES

Medical device	No		
Marquage	CE0123	Barcode	3661040217963
Guarantee (month):	12	Customs Code	901890
Sterile	No		

INFORMATIONS

Main unit

Packaged unit

Width (cm) :	Width packaging (cm):
Length (cm) :	Length packaging (cm) :
Height (cm) :	Height packaging (cm) :
Weight (kg) :	Weight packaging (kg) :

SPECIFIC STORAGE CONDITION

Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:47 - This article is available on our website : www.gsh-med.fr with the following reference DS99325_N

All characteristics are given for information only and may be changed