

SUCCION UNIT NEW ASKIR 118 BASIC-WM, 1 LITER FLOVAC DISPOSABLE LINER

Réf. BCASPI178

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
II
a



SUCCION UNIT NEW ASKIR 118 BASIC-WM, 1 LITER FLOVAC DISPOSABLE LINER



Réf. BCASPI178

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe IIa | |
| Marquage | CE0123 | Barcode | 8054619104712 |
| Guarantee (month): | 24 | Customs Code | 90189084 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 19,0 | Width packaging (cm): | |
| Length (cm) : | 35,0 | Length packaging (cm) : | |
| Height (cm) : | 15,0 | Height packaging (cm) : | |
| Weight (kg) : | 2,500 | Weight packaging (kg) : | |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 03:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI178