

## NC3 MONITOR NIBP+SPO2 COMEN

COMEN

Réf. BCMONI01



Batter



Secteur

CE  
1639

DISPOSITIF MEDICAL CLASSÉ II  
MEDICAL DEVICE



# NC3 MONITOR NIBP+SPO2 COMEN

COMEN

Réf. BCMONI01

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE1639    | Barcode                 | 3661040219103 |
| Guarantee (month):         | 36        | Customs Code            | 9018199560    |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 21,0      | Width packaging (cm):   | 25,5          |
| Length (cm) :              | 13,0      | Length packaging (cm) : | 31,5          |
| Height (cm) :              | 13,0      | Height packaging (cm) : | 26,5          |
| Weight (kg) :              | 1,250     | Weight packaging (kg) : | 2,950         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -20       | Temperature max         | 60            |
| Humidity min               |           | Humidity max            | 93            |



07/11/2022 12:34 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCMONI01

All characteristics are given for information only and may be changed