

NC5 MONITOR NIBP+SPO2 COMEN+TEMP+ECG+PRINT

COMEN

Réf. BCMONI05



Batterie



Secteur

CE
1639

DISPOSITIF MEDICAL
MEDICAL DEVICE CLASS II



NC5 MONITOR NIBP+SPO2 COMEN+TEMP+ECG+PRINT

COMEN

Réf. **BCMONI05**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE1639	Barcode	3661040219097
Guarantee (month):	36	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	16,5	Width packaging (cm):	29,5
Length (cm) :	16,5	Length packaging (cm) :	36,5
Height (cm) :	25,0	Height packaging (cm) :	26,0
Weight (kg) :	2,300	Weight packaging (kg) :	3,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min		Humidity max	85



17/05/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference BCMONI05

All characteristics are given for information only and may be changed