

NC3 MONITOR NIBP+SPO2 COMEN+TEMP

COMEN

Réf. BCMONI02




Batteries



Secour

CE
1639

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE II



NC3 MONITOR NIBP+SPO2 COMEN+TEMP

COMEN

Réf. BCMONI02

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE1639	Barcode	3661040219066
Guarantee (month):	36	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,5	Width packaging (cm):	25,5
Length (cm) :	13,0	Length packaging (cm) :	31,5
Height (cm) :	22,0	Height packaging (cm) :	26,5
Weight (kg) :	1,250	Weight packaging (kg) :	2,950
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min		Humidity max	85



07/11/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference BCMONI02