

SPO2 SILICON WRAP SENSOR FOR COMEN NC3/NC5/C34

Réf. DBSWRAP01



SPO2 SILICON WRAP SENSOR FOR COMEN NC3/NC5/C34

Réf. **DBSWRAP01**

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0537	Barcode	3661040219455
Guarantee (month):	12	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	13,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	50
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:36 - This article is available on our website : www.gsh-med.fr with the following reference DBSWRAP01

All characteristics are given for information only and may be changed

2/2