

## REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

Réf. BCTHER01



# REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

Réf. **BCTHER01**

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE1639    | Barcode                 | 3661040219349 |
| Guarantee (month):         | 6         | Customs Code            | 9018193010    |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              |           | Weight packaging (kg) : |               |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -25       | Temperature max         | 55            |
| Humidity min               |           | Humidity max            | 85            |



31/03/2022 03:34 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCTHER01

All characteristics are given for information only and may be changed