

## NC5 MONITOR NIBP + SPO2 COMEN + TEMP + PRINT

COMEN

Réf. BCMONI09



Batterie



Secteur

CE  
1639

DISPOSITIF MEDICAL  
MEDICAL DEVICE CLASS II



# NC5 MONITOR NIBP + SPO2 COMEN + TEMP + PRINT

COMEN

Réf. **BCMONI09**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE1639	Barcode	3661040220017
Guarantee (month):	36	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	16,5	Width packaging (cm):	29,5
Length (cm) :	16,5	Length packaging (cm) :	36,5
Height (cm) :	25,0	Height packaging (cm) :	26,0
Weight (kg) :	2,300	Weight packaging (kg) :	3,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	50
Humidity min		Humidity max	85



17/05/2022 12:57 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCMONI09

All characteristics are given for information only and may be changed