

## HIGH VISIBILITY EMT VEST ELITE BAGS S/M - E-VEST'S - RED EB02.055

 ELITE BAGS

Réf. SEGILE01

GARANTIE  
**2**  
ANS  
YEARS  
GUARANTEE



## HIGH VISIBILITY EMT VEST ELITE BAGS S/M - E-VEST'S - RED EB02.055



Réf. **SEGILE01**

FEATURES			
Medical device	No		
Marquage	Barcode		
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	59,0	Width packaging (cm):	3,0
Length (cm) :	54,0	Length packaging (cm) :	59,0
Height (cm) :	3,5	Height packaging (cm) :	54,0
Weight (kg) :	0,970	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/06/2022 12:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SEGILE01

All characteristics are given for information only and may be changed