

## CHARGING BASE STATION FOR THERMOSCAN PRO 6000 TYMPANIC THERMOMETER

Réf. DWPRO6K07

WelchAllyn®



# CHARGING BASE STATION FOR THERMOSCAN PRO 6000 TYMPANIC THERMOMETER



Réf. **DWPRO6K07**

FEATURES			
Medical device	No		
Marquage	CE	Barcode	00732094325324
Guarantee (month):	36	Customs Code	850410
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	16,0
Length (cm) :		Length packaging (cm) :	25,0
Height (cm) :		Height packaging (cm) :	9,0
Weight (kg) :		Weight packaging (kg) :	0,700
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:56 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DWPRO6K07