

## CHAIR SCALE EKMCC01, CLASS III VERIFICATION

Réf. EKMCC01\_CIII

**KERN**



CE  
0297



# CHAIR SCALE EKMCC01, CLASS III VERIFICATION



## Réf. EKMCC01\_CIII

Mobile version with locking brakes on all wheels

- With its four wheels, this chair scale provides maximum mobility in bringing the scale to the patient. This ensures more efficient use of time for the clinic staff and greater safety for the patients, who can be weighed in their familiar environment.
- This chair scale is the ideal measuring instrument for retirement homes, rehabilitation centers and clinics with lots of patients, and can handle weighing overweight persons weighing up to 250 kg
- For fragile and weak patients, the comfortable chair design offers safe, secure seating during weighing
- Foldable armrests and footrests make transfers into the chair easier. Ideal for overweight patients or barrier-free use, e. g. for transfers from a bed to the scale
- Ergonomically positioned carrying handles
- Automatic Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically «frozen» for 10 sec. This means that you have sufficient time to attend to the patient first and then take the weight reading in peace
- BMI function to determine normal weight / surplus weight
- High resolution readout: Read-out [d] for 5 sec increased by one decimal place at the touch of a key
- Practical mains adapter pouch to store the mains adapter which is supplied as standard. Fitting option using two Velcro fasteners.
- Large backlit LCD display, digit height 25 mm

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement
Marquage	CE0297		Classe III
Guarantee (month):	36	Barcode	3661040055060
Sterile	No	Customs Code	
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :	99,0	Width packaging (cm):	
Length (cm) :	62,5	Length packaging (cm) :	
Height (cm) :	98,5	Height packaging (cm) :	
Weight (kg) :	24,300	Weight packaging (kg) :	37,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:04 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMCC01\_CIII

All characteristics are given for information only and may be changed

## CHAIR SCALE EKMCC01, CLASS III VERIFICATION



Réf. EKMCC01\_CIII

### COMPOSITION

Code	Descriptive	Quantity
EKMCC01	Chair scale Kern MCC - 250 Kgs	1
EKCERTI06	EU-Verification (CE) electr.balances 50 kg - 350 kg	1
EKTRANSPORT02	Freight Cost Big Products	1