

THERMAL RECEIPT ROLLS FOR PRINTER KERN YKB-01N (10 PIECES)

KERN

Réf. EKMWA08



THERMAL RECEIPT ROLLS FOR PRINTER KERN YKB-01N (10 PIECES)



Réf. **EKMWA08**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040173719
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	20,0
Length (cm) :		Length packaging (cm) :	8,0
Height (cm) :		Height packaging (cm) :	11,0
Weight (kg) :		Weight packaging (kg) :	1,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:34 - This article is available on our website : www.gsh-med.fr with the following reference EKMWA08