

KERN MPE SCALE WITH HEIGHT 250 KG WITH CLASS III VERIFICATION

Réf. EKMPE02-CIII

KERN



CE
0297



KERN MPE SCALE WITH HEIGHT 250 KG WITH CLASS III VERIFICATION



Réf. EKMPE02-CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297	Barcode	3661040063249
Guarantee (month):	36	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	45,0	Width packaging (cm):	
Length (cm) :	36,5	Length packaging (cm) :	
Height (cm) :	103,0	Height packaging (cm) :	
Weight (kg) :	13,000	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:05 - This article is available on our website : www.gsh-med.fr with the following reference EKMPE02-CIII

KERN MPE SCALE WITH HEIGHT 250 KG WITH CLASS III VERIFICATION



Réf. EKMPE02-CIII

COMPOSITION

Code	Descriptive	Quantity
EKMPE02	Kern MPE Scale with height - 250 kg	1
EKCERTI06	EU-Verification (CE) electr.balances 50 kg - 350 kg	1