

## S/S THERMOMETER JAR H.100 X 30 (DIAM) MM

Réf. AI11611



## S/S THERMOMETER JAR H.100 X 30 (DIAM) MM



Réf. AI11611

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019291
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,0	Width packaging (cm):	8,0
Length (cm) :	10,0	Length packaging (cm) :	11,5
Height (cm) :	10,0	Height packaging (cm) :	8,0
Weight (kg) :	0,088	Weight packaging (kg) :	0,088
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:44 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI11611

All characteristics are given for information only and may be changed