

## INTERFACE CABLE RS 232 FOR KERN MWS, MXS, MTS, MPS

**KERN**

Réf. EKMP508



## INTERFACE CABLE RS 232 FOR KERN MWS, MXS, MTS, MPS



Réf. **EKMPS08**

FEATURES			
Medical device	No		
Marquage		Barcode	4045761096950
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMPS08

All characteristics are given for information only and may be changed