

S/S KOCHER FORCEPS JAR 25 (DIAM)X H. 140 MM

Réf. AI11615



S/S KOCHER FORCEPS JAR 25 (DIAM)X H. 140 MM



Réf. AI11615

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019307
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,5	Width packaging (cm):	9,0
Length (cm) :	14,0	Length packaging (cm) :	17,0
Height (cm) :	14,0	Height packaging (cm) :	9,5
Weight (kg) :	0,121	Weight packaging (kg) :	0,126
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:45 - This article is available on our website : www.gsh-med.fr with the following reference AI11615

All characteristics are given for information only and may be changed