


SURGERY SET WITH ALUMINIUM BOX 18X9X3 CM

Réf. IA05601



SURGERY SET WITH ALUMINIUM BOX 18X9X3 CM



Réf. **IA05601**

Consisting of:

- IC09411 Straight iridectomy scissors, 11.5 cm
- IC11414 Straight round-ended scissors, 14 cm
- ID05403 no.3 scalpel handle
- IP02014 Straight Kocher toothed forceps, 14 cm
- IP21414 Pean Murphy forceps, 14 cm
- IP26014 toothed dissection forceps, 14 cm
- IP26014 Dissection forceps without teeth, 14 cm
- IP26410 Splinter forceps, 10 cm
- IP29414 Mayo-Hegar needle holder, 14 cm
- IS09414 Cotton holding stylet, 14 cm

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040006314
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,0	Width packaging (cm):	
Length (cm) :	18,0	Length packaging (cm) :	
Height (cm) :	3,0	Height packaging (cm) :	
Weight (kg) :	0,450	Weight packaging (kg) :	0,450
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:09 - This article is available on our website : www.gsh-med.fr with the following reference IA05601

All characteristics are given for information only and may be changed

SURGERY SET WITH ALUMINIUM BOX 18X9X3 CM



Réf. IA05601

COMPOSITION

Code	Descriptive	Quantity
AAT1893	Aluminium Box Grey 18x9x3 cm	1
IA05602	Minor surgery set for IA05601	1