

INFANT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT2(13 TO 20 CM)

Réf. DS99506



CE
0483



INFANT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT2(13 TO 20 CM)



Réf. **DS99506**

Infant cuff (13 to 20 cm)
Washable nylon cuff without pocket with Velcro fastening

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0483	Barcode	3661040075259
Guarantee (month):	12	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	11,2	Length packaging (cm) :	
Height (cm) :	2,3	Height packaging (cm) :	
Weight (kg) :	0,052	Weight packaging (kg) :	0,052
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	10	Humidity max	85



30/03/2022 14:59 - This article is available on our website : www.gsh-med.fr with the following reference DS99506

All characteristics are given for information only and may be changed