

S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM

Réf. AI14830



S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM



Réf. AI14830

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019352
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	30,0	Width packaging (cm):	55,0
Length (cm) :	48,0	Length packaging (cm) :	30,0
Height (cm) :		Height packaging (cm) :	8,5
Weight (kg) :	1,070	Weight packaging (kg) :	1,080
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:45 - This article is available on our website : www.gsh-med.fr with the following reference AI14830

All characteristics are given for information only and may be changed