

LID FOR BEDPAN (STAINLESS STEEL) 48X30CM

Réf. AI14831



LID FOR BEDPAN (STAINLESS STEEL) 48X30CM



Réf. AI14831

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | CE | Barcode | 3661040032917 |
| Guarantee (month): | 12 | Customs Code | 90189090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 30,0 | Width packaging (cm): | 54,0 |
| Length (cm) : | 48,0 | Length packaging (cm) : | 39,0 |
| Height (cm) : | | Height packaging (cm) : | 2,0 |
| Weight (kg) : | 0,416 | Weight packaging (kg) : | 0,416 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 14:45 - This article is available on our website : www.gsh-med.fr with the following reference AI14831

All characteristics are given for information only and may be changed