


## S/S BOWL WITH SPOUT DIAM 40 MM - 20 ML

Réf. AI15040



## S/S BOWL WITH SPOUT DIAM 40 MM - 20 ML



Réf. AI15040

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019369
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,0	Width packaging (cm):	5,0
Length (cm) :		Length packaging (cm) :	5,0
Height (cm) :		Height packaging (cm) :	2,5
Weight (kg) :	0,025	Weight packaging (kg) :	0,025
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/07/2022 12:46 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI15040

All characteristics are given for information only and may be changed