


S/S BOWL WITH SPOUT DIAM 240 MM - 3LITERS

Réf. AI15240



S/S BOWL WITH SPOUT DIAM 240 MM - 3LITERS



Réf. AI15240

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019444
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	24,0	Width packaging (cm):	24,7
Length (cm) :	24,0	Length packaging (cm) :	24,7
Height (cm) :	11,6	Height packaging (cm) :	11,6
Weight (kg) :	0,241	Weight packaging (kg) :	0,241
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:46 - This article is available on our website : www.gsh-med.fr with the following reference AI15240

All characteristics are given for information only and may be changed