

S/S BOX 18 X 8 X 4 CM

Réf. AI18084



S/S BOX 18 X 8 X 4 CM



Réf. AI18084

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019482
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	19,0
Length (cm) :	18,0	Length packaging (cm) :	21,5
Height (cm) :	4,0	Height packaging (cm) :	10,0
Weight (kg) :	0,256	Weight packaging (kg) :	0,265
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:47 - This article is available on our website : www.gsh-med.fr with the following reference AI18084

All characteristics are given for information only and may be changed