

S/S BOX 20 X 10 X 3 CM - TOP/BOTTOM PERFORATED

Réf. AI20103P



 Holtex



S/S BOX 20 X 10 X 3 CM - TOP/BOTTOM PERFORATED



Réf. AI20103P

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040067735
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	13,6	Width packaging (cm):	13,6
Length (cm) :	28,0	Length packaging (cm) :	28,0
Height (cm) :	5,5	Height packaging (cm) :	5,5
Weight (kg) :	0,325	Weight packaging (kg) :	0,325
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:47 - This article is available on our website : www.gsh-med.fr with the following reference AI20103P

All characteristics are given for information only and may be changed