

S/S BOX 20 X 10 X 5 CM

Réf. AI20105



S/S BOX 20 X 10 X 5 CM



Réf. AI20105

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019505
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,6	Width packaging (cm):	11,6
Length (cm) :	20,5	Length packaging (cm) :	21,6
Height (cm) :	5,1	Height packaging (cm) :	6,8
Weight (kg) :	0,397	Weight packaging (kg) :	0,397
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:47 - This article is available on our website : www.gsh-med.fr with the following reference AI20105

All characteristics are given for information only and may be changed