

S/S BOX PERFORATED 25 X 12 X 6 CM PERFORATED TOP / BOTTOM

Réf. AI25126P



S/S BOX PERFORATED 25 X 12 X 6 CM PERFORATED TOP / BOTTOM

Réf. AI25126P

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040067773
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,9	Width packaging (cm):	14,5
Length (cm) :	25,7	Length packaging (cm) :	26,5
Height (cm) :	6,6	Height packaging (cm) :	8,2
Weight (kg) :	0,613	Weight packaging (kg) :	0,670
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference AI25126P